



Frequently Asked Questions

Q. What is the Transformational Beauty Skin Rejuvenation System?

A. An at-home wrinkle-reducing system clinically proven to reduce fine lines and wrinkles around the eyes, mouth and brow.

Q. What are the advantages of using the Transformational Beauty System?

A. The Transformational Beauty System goes deeper than typical creams and serums to dramatically reduce the appearance of fine lines and wrinkles — in the privacy and comfort of you own home. No painful needles. No costly appointments. It is simple and easy to use.

Q. How do I know if the Transformational Beauty System is right for me?

A. If you have fine lines and wrinkles around your eyes, mouth or forehead, then the System is for you. It is appropriate for adult women and men of all ages and skin types.

Q. How do I use the Transformational Beauty System?

A. In 4 easy steps:

- Step 1 - Use the Gentle Exfoliating Cream to remove dead skin and prepare skin for optimal treatment.
- Step 2 – Clean skin with Power Infusion Skin Prep
Clean the face and remove any makeup, moisturizers, or other skin products. Wait for skin to dry. IT'S VERY IMPORTANT TO MAKE SURE SKIN IS CLEAN AND DRY.
- Step 3 – Treat with Power Infusion Eye, Brow or Lip Patch.
Remove patches from foil. Place on countertop and attach both buds. Apply Eye, Brow or Lip patch and wear for 30 minutes. Treatment will start and stop automatically.
- Step 4 - Maintain with Youth Serum HA
Remove patch slowly while pressing lightly on skin to prevent pulling. Apply Youth Serum HA to seal in the treatment for long lasting results.

Q. Will I feel the Transformational Beauty System working?

A. Yes. You may feel a slight tingle that moves from one patch to the other.



Q. How often do I need to use the Transformational Beauty System?

A. For optimal results, use the Transformational Beauty System twice a week during the first week, then maintain your results with weekly treatments (once a week).

Q. How long does each treatment session last?

A. Each treatment lasts just 30 minutes and can be completed at home, on your own time.

Q. How quickly should I see results?

A. Many individuals experience a significant reduction in the appearance of wrinkles with just one treatment. Speed of results varies depending on severity of wrinkles, and individual skin factors. Allow 4 weeks to see the accumulative results and benefits of consistent use.

Q. Can I apply more than one set of patches on the same day?

A. No. Do not apply patch treatments more than once every three days. Use only as directed.

Q. Can I reuse the patches?

A. No. Patches are only for single use.

Q. Should I be doing anything else during treatments?

A. To help prevent accidental damage to the Power Pod, avoid too much activity during treatments. For best results, facial muscles should be in a relaxed state.

Q. Can I apply other products on the patches?

A. No. The ingredients are already infused into each patch. Never apply any other products on to the patches.

Q. What if the patches do not stick properly?

A. Make sure the skin around the treatment area is completely clean and dry before applying each patch treatment. The patches will adhere to the skin as long as the treatment area is thoroughly cleansed and dry. Recommended: Alcohol-free Power Perfect Skin Prep effectively and gently cleanses, plus improves connectivity for micro-charged infusion.



Q. How do I dispose of the patches?

A. The patches may be disposed of as household waste. **Tip:** Once your treatment is finished, remove the patches from your face, turn them gel-side up and stick them together. Then remove the Power Pod Connector Buds.

Q. What if the Power Pod cords are too snug around the neck?

A. If needed, slide down the small connector on the cords until it's about an inch above the Power Pod. Adjust the fit by gently separating the cords from each other. Be very careful not to pull too hard; they should separate easily.

Q. Is it important to tuck the Power Pod cords behind my ears?

A. For optimal results, tuck the Power Pod cords behind your ears. This will help prevent the patches from pulling on the skin.

Q. I snapped the Power Pod buds onto each patch and applied them but the pod isn't working. How do I reset it?

A. On the rare occasion that your Power Pod locks up, simply unsnap the Power Pod buds from the patches (you can keep the patches on your face) and touch the two buds together until a green light appears on the Power Pod. Reapply the buds onto each patch and resume your treatment.

Q. How does the Transformational Beauty System work?

A. The Transformational Beauty System uses patented infusion technology to deliver Hyaluronic Acid — the wrinkle-filler injected by dermatologist — and peptides by infusing them into the top layers of your skin, unlike traditional creams or serums that only sit on top of the surface of the skin.

Q. Was the Transformational Beauty System tested for safety?

A. Yes, the Transformational Beauty System has been safety-tested with hundreds of clinical studies participants.

Q. If I am pregnant or breastfeeding can I use the System?

A. We recommend you speak with your physician.

Q. I have metal implants in my mouth, can I use the System?

A. Yes, however you might experience a slight metallic taste which is not harmful.

Q. I have epilepsy, can I use the System?

A. No.



Q. I have a pacemaker, can I use the system?

A. No.

Q. Can the Transformational Beauty System be purchased or shipped internationally?

A. The Transformational Beauty System is only available for sale and shipment within the United States, Alaska, Hawaii and Puerto Rico.

Q. What if I buy the Transformational Beauty System and am not satisfied, Can I get a refund?

A. Yes. The Transformational Beauty System comes with a 30-day Money Back Guarantee. If you decide you do not want to keep the Transformational Beauty System for any reason within 30 days of your purchase, view our Return Policy and follow the instructions.

Q. I've lost the Transformational Beauty System brochure/instructions. How do I get another?

A. Please call Customer Service at 1.855.777.0245, Monday—Friday, 8 a.m.–5 p.m. (PST)

Q. I've misplaced the Transformational Beauty System Power Pod. How do I get another?

A. Please contact Customer Service at 1.855.777.0245, Monday—Friday, 8 a.m.–5 p.m. (PST).

PRECAUTIONS:

Keep the Power Pod and patches out of reach of children. Use only as directed. Check with your physician before using if you are pregnant, have a pacemaker or similar device, are epileptic, wear metal braces on your teeth, or have a metal implant. Do not use over an area with broken capillaries or over an inflamed or open wound. Do not disassemble the Power Pod. The rechargeable battery does not require access. The pod contains no reusable products. Do not recharge, puncture, modify, dispose in fire or trash incinerator, or heat above 212 °F (100 °C). Do not submerge in water.